


Hello, **Wireless**



ArthroFree[®]

Enhanced Imaging for the
Wireless Age of Surgery

With **ArthroFree 2**, experience a crystal clear, crisp image with **Meridiem[®] image enhancement** coupled with the freedom to move between ports without being tied to the surgical tower. As your camera hand guides your vision, complete your procedures with unparalleled focus and efficiency. Have a full day in the OR without cables dragging down your hands with every move.




Experience the **Freedom of Wireless**

Focus

The ArthroFree 2 System:


- ✓ Eliminates camera-cable management & interruptions
- ✓ Enables you to fix your eyes on the screen—and keep them there
- ✓ Delivers a crystal clear, crisp surgical image



Efficiency

The ArthroFree 2 System:

- ✓ Reduces set-up and take-down time
- ✓ Integrates seamlessly with existing equipment
- ✓ Eliminates slip, trip, fall and burn hazards in the OR



ArthroFree[®]

S4K Wireless Operating Room Kit



- 1 32" 4K Surgical Monitor**

The 4K UHD surgical monitor faithfully reproduces 4K endoscope images in Ultra High Definition. With a 178° viewing angle and technology to reduce unwanted reflections and glare, experience a highly visible surgical image from anywhere in the operating field.
- 2 MVR 4K Recorder**

Next generation DICOM and AI ready recorder with embedded neural processor. 2TB 2.5"HDD for faster capture and processing. 7" touch screen for ease of navigation.
- 3 Meridiem[®] Image Enhancer**

Bringing image clarity and sharpness control and adjustment to the Operating Room, the system allows surgeons to optimize the visual display for each unique surgical situation, ensuring precision and confidence in every procedure.
- 4 ArthroFree Wireless Camera Head**

1080p HD camera, featuring the on-board reimagined Meridiem[®] Light Engine, delivers 3x more light for enhanced visibility. Paired with optimized image processing firmware and an improved white balance algorithm, it ensures a crisp, clear surgical image.
- 5 Wireless Receiver**

Wireless transfer from the camera to the receiver using encrypted ultra-wide band (UWB) technology that does not rely on customer IT systems. The system delivers smooth real-time video ensuring optimal clarity and precision for surgical visualization.
- 6 Charger and Batteries**

Guaranteed to last, the battery delivers 60 min of power even after 300 charges. A 5-bay charger ensures every case is ready to go the distance, keeping your workflow uninterrupted.

SCHEDULE A DEMO TODAY

www.lazurite.co | 1.833.214.2324

sales@lazurite.co